


Standard Patient Transfer Slings




- sling shapes -



















- measuring -














- training notes -

Fastfit Sling

Optional Extras:

Head Support, Hip Tapes, **Chest/Waist Support**

The Fastfit sling is a good general use sling.

It is a divided leg style sling and has padded legs for user comfort. The legs can be 'hammocked' together if needed. It has a large enough commode aperture to facilitate toileting if needed.

Fastfit Deluxe Sling

Optional Extras:

Head Support, Hip Tapes, **Chest/Waist Support**

A good general use sling. It picks up from the outside of the thighs and supports the entire thigh area down to the knee. Ideal for people with painful hips. The leg sections can also be 'hammocked' if necessary to avoid coming between the

A more reclined position may be better for the hammocked pick up if used that way. Can be used for long limb amputations (knee or below).

Universal Deluxe Sling

Optional Extras:

Head Support, Hip Tapes, **Chest/Waist Support**

A very good style of sling for users with hip problems or who require additional support.

It has padded legs and supports each thigh fully.

The optional hip tapes provide a more even weight distribution and take pressure from behind the knee.

Longseat Sling

Optional Extras:

Head Support, Hip Tapes, **Chest/Waist Support**

A very basic design of sling which has to be applied on the bed.

Good for some amputees and users who do not want the indignity of a sling between their legs.

Can have a commode aperture if required.

The padded roll at the bottom of the sling will provide a more comfortable transfer for users.

Dress/Toileting Sling

Optional Extras:

Head Support

This sling is used for transfers to and from a toilet or commode.

It has a large aperture and therefore the user must have good upper body tone in order to use this sling. A padded underarm roll provides additional comfort during the transfer. Also available without padding. It has a Velcro waist strap with a clunk-click buckle for additional safety.

Optional Extras:

Chest/Waist Support, **Loop Fixings, Dual Fixings** (clips & loops)

A very good cocooning style sling.


High Easy Sling


Can only be used for stud fixing hoists. It has an integral head support as the user is hoisted in a slightly reclined position. Especially good for smaller children due to the cocooning nature of the sling.


Only sling available for stud fixing type hoists although a toileting version is available if needed.


Does not suit all types of disability but you are limited if the hoist is stud fixing.


- optional extras -














CPD TRAINING COURSES AVAILABLE ONLINE OR ON-SITE

Global Training by Silvalea have developed a series of specific training courses designed to deliver high quality, information led training to professionals. The courses feature a structured and methodical approach to learning about the chosen subject, whether it be safety regulations or measuring an end user. Based on the structure and value of our training courses, the CPD approved that Silvalea adheres to the full conformity of the CPD guidelines.

www.silvalea.com